

QuantiVirus™ Saliva Collection Kit

For COVID-19 testing



Easy sampling for both patients and the physician



Samples can be stable for 7 days in sample preservation solution



Samples can be directly use for RNA isolation

DiaCarta provides limited supply of high-quality saliva collection kits that meet COVID-19 testing needs.

The QuantiVirus™ Saliva Collection Kit is a sample collection kit of patient saliva for SARS-CoV-2 and other virus detection. The kit uses a simple procedure for saliva collection, provides comfort to the patient, and maintains sample stability up to 7 days at room temperature before use for testing.

Product Specification and Ordering Information

Product Name	QuantiVirus™ Saliva Collection Kit
Packing Size	1 sample per kit
Catalog Number	DC-11-00021
Kit Components	The QuantiVirus™ Saliva Collection Kit contains a collection funnel, tube and a storage tube containing sample preservation solution. The storage tube is sealed in a biohazard bag for shipping or processing.
Storage and Shipping Instruction	<p>Before use, the kit can be shipped and stored at room temperature for 24 months</p> <p>After the kit is used, the sample can be shipped stored at room temperature for up to 7 days before processing. Samples must be processed within 7 days of saliva collection.</p>